d.Flow

Flowmeter


Am Parir 27
52379 Langerwehe, Germany
Phone: +49 2423 / 40493-5
Mail: info@ziemer-chromatographie.de
www.ziemer-chromatographie.de

Flowmeter d.Flow


The flowmeter is especially designed for checking the accuracy and fluctuation of HPLC or dosing pumps. The device is calibrated and traceable according to DKD or NIST. The flowmeter is supplied with a power adapter, tripod, test certificate and transport case.

Figure 1: Operating schematic

Operating principle

Figure 2: Measuring cycle

The flowmeter d.Flow determines the actual flow rate of a liquid flow. For this purpose, the time to fill a measuring tube with constant volume is determined. The display shows the measured value and the filling level. The measured values could be logged or transmitted directly via RS232.

The measurement is volumetric by opto-electronic detection of the liquid level in a rising pipe. The time that passes between reaching a lower and an upper level is being measured and the flow is being calculated. After the upper level has been reached, the rising pipe will be drained via a magnetic valve and the measurement can start again.

Density and temperature of the medium do not matter due to the volumetric measuring principle. Due to the low surface tension, water, ethanol or methanol are recommended.

Advanced features

Detailed measurement protocol

Printing function

Data logger

Data transfer to PC

```
*** Calibration Report ***  
Date/Time    22-09-2021 08:22:38  
  
Tester Company  
             DURATEC Analysentechnik  
Tester Name   Thomas Nied  
Flowmeter Type  DURATEC dFlow  
Flowmeter Serial 1840052140  
  
Pump Manufacturer  DURATEC  
Pump Model         dDrive C30  
Pump Serial        4711  
Pump Flow [ml/min] 2.00  
  
--Time--           --Flow--  
08:22:40           1.9817 ml/min  
08:22:53           1.9898 ml/min  
08:23:06           1.9894 ml/min  
08:23:19           1.9830 ml/min  
08:23:32           1.9929 ml/min  
08:23:45           1.9910 ml/min  
08:23:58           1.9892 ml/min  
08:24:11           1.9898 ml/min  
08:24:24           1.9866 ml/min  
08:24:36           1.9871 ml/min  
  
-----  
Signature
```

Calibration certificate / recalibration

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vFlow Kalibrierfaktor 13542
vFlow Calibration Factor

The flowmeter includes a calibration certificate. The flowmeter is checked at 3 specified flow rates. Each setpoint is measured 10 times. The test equipment used for calibration is traceable to DKD. Optionally, it is possible to use up to 6 preset points. In order to be able to use the flowmeter as a test equipment without any time restrictions, an annual recalibration is recommended.

Specification

Packing list	1x Flowmeter d.Flow with tripod, power supply, waste cup, manual, certificate (with 3-point calibration), transport case
Fluid path	PEEK, PTFE, FEP, EPDM
Flow range	0.05 ml/min - 20 ml/min
Measuring principle	volumetric
Measuring volume	approx. 240 µl
Accuracy	better 1% (measuring range 0.1 - 10.0 ml/min)*
Reproducibility	better 1% (measuring range 0.1 - 10.0 ml/min)*
Resolution	0.0001 ml/min
Display	3.5" touch display (LCD-TFT)
Output	RS232 interface
Power supply	24 V/DC
Dimensions	Flowmeter: W: 105 mm / H: 170 mm / D: 40 mm (without tripod) Case: W: 420 mm / H: 360 mm / D: 135 mm
Weight	approx. 3.1 kg (incl. case and accessories)

*based on water as reference medium

Ordering Information

Part no.	Description
792625	Flowmeter d.Flow with tripod, power supply, waste cup, manual, certificate (with 3-point calibration), transport case
792615-28	Thermoprinter for flowmeter with power supply and RS232 connection cable
792615-CAL	3-point recalibration incl. certificate
792615-CAL-4	4-point recalibration incl. certificate
792615-CAL-5	5-point recalibration incl. certificate
792615-CAL-6	6-point recalibration incl. certificate

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DURATEC Analysetechnik GmbH