

HPLC Multo - Phases

Phase	Bonded Phase	Pore Size	Available Particle Size	Comments
CS-ODS 100	C18	100 Å	5µ	Cheap Standard Phase
MultoHigh-U	C18, C18 AQ	120 Å	2,1µ, 1,7µ (C18)	UHPLC Phase
MultoHigh	C18, C8, Si, Phenyl, Diol, Cyano, Amino	100 Å	3µ, 5µ, 10µ	
MultoHigh-Bio	C18, C8	200 Å	5µ, 10µ	Bio-Separation Phase
	C4	300 Å	5µ, (10µ)	Bio-Separation Phase
Multospher	C18, C8, Si	100 Å	3µ, 5µ, 10µ	Standard Phase
	C18, C8, Si	120 Å	3µ, 5µ, 10µ	Standard Phase
	C18-AQ, C18-HP	120 Å	3µ, 5µ	High Purity Phase
	C18-FBS	100 Å	5µ	Polymeric Bonded for Acids and Bases
	PAH-III	300 Å	5µ	for Poly Aromatic Hydrocarbons (PAH)
	PSM	100 Å	5µ	for Pesticides
MultoKrom	C18, C8, C4, Si	100 Å	3µ, 5µ, 10µ	Silica Specification analog Kromasil
Multohyp	C18, BDS-C18	120 Å	3µ, 5µ	Silica Specification analog Hypersil
MultoWorb	C18, C8, C6, C1, Si, Phenyl, CN, SCX, SAX, NH2	80 Å	3µ, 5µ	Silica Specification analog Spherisorb



High Quality stainless steel columns.
Refill-Service available.



High Quality stainless steel cartridges as guard columns.
(Compatible with CS-type columns from 2 – 4.6 mm)